

ADA Specifications for Pool Lifts

Specification Description	Technical Specification	RMT Aquatics Lifts
Pool Lift Location	Maximum Water Depth of 48 inches	All lifts will operate in this environment
Clear Deck Space	36 inch wide by 48 inch deep transfer area that starts 12 inches behind the seat	All lifts will operate in this environment. PAL Lifts equipped with 240 degree turning radius. Splash Lifts have 359 degree turning radius
Seat Height	The seat must be able to stop between 16 inches and 19 inches above the pool deck.	All Lifts are capable of stopping at specified distance
Seat Width	Seat width minimum of 16 inches	Our seats are 18.5 inches wide
Foot and Arm Rests	Footrests are required. Armrests are optional, but must be foldable or removable when transferring.	Footrests are furnished with all lifts. Detachable, foldable armrests are available as an optional accessory for all lifts.
Operation	Lift must be capable of unassisted operation. A person must be able to call the lift when it is in either the deck or water position.	All lifts come with a completely waterproof remote control that is attached to the unit. This control will operate from either the deck or water location.
Submerged Depth	Lift must submerge to 18 inches below water surface minimum.	All lifts will submerge in excess of 18 inches when installed properly.
Lifting Capacity	Minimum of 300 lbs	All lifts will lift at least 300 lbs. Splash HD lift will transfer users up to 400lbs

If you have any questions, or would like more information, please call us at (800) 577-4424.